RISK GROUP 1: NOTICE IR EMITTED FROM THIS PRODUCT.
USE APPROPRIATE SHIELDING OR EYE PROTECTION.

©2026 STYLIDEAS Ltd. All rights reserved.
STYLIDEAS owns trademarks, registered trademarks, patents, patents pending, design
registrations and protection of the copyright of its products. www.stylideas.com/ip.
Distributed by STYLIDEAS Ltd.
UK: Amshold House, IG10 2RW, UK. +44(0)20 3225 5559 teamtom@stylideas.com.
EU RP: OBELIS S.A. Bd. Général Wahnis 53, B-1030 Brussels, Belgium


RBP01A
Made in China.

STYLPRO™

Red Light Full Body Panel

Panneau de lumière rouge pour tout le corps
Rötlicht-Panel für den ganzen Körper
Rood licht-plaat voor het hele lichaam
Pannello corpo con luce rossa
Panel do nasświetlania całego ciała światłem czerwonym
Alfombrilla de terapia de luz roja para todo el cuerpo
Rödljuspanel för hela kroppen
Pannello corpo con luce rossa
Panel do nasświetlania całego ciała światłem czerwonym
Alfombrilla de terapia de luz roja para todo el cuerpo
Rödljuspanel för hela kroppen
Pannello corpo con luce rossa
Panel do nasświetlania całego ciała światłem czerwonym
Alfombrilla de terapia de luz roja para todo el cuerpo
Rödljuspanel för hela kroppen


Quick Guide


How To Use

***FOR YOUR SAFETY AND CONTINUED ENJOYMENT OF THIS PRODUCT, ALWAYS READ INSTRUCTION MANUAL CAREFULLY BEFORE USING THE APPLIANCE.**

1. Plug in the device.

2. Always wear the provided eye protection glasses.

3. Position or wrap the panel around the desired area.


4. Long press the power on/ off button to start the session.

5. Relax. Panel will auto off after 20 minutes.

Using the Panel


Hold the On/Off button to turn on the panel. Four display lights will appear near the On/Off button, each representing 5 minutes. After every 5 minutes of use a light will turn off. The panel will automatically turn off after 20 minutes.

Warning

Do not fold or sit on the panel. Always wear eye protection. Closing your eyes will not protect them from invisible NIR light.

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Strap Guide

How To Wear

*THE STRAP GUIDE BELOW OUTLINES THE RECOMMENDED USES AND APPROPRIATE STRAPS TO USE FOR EACH WEARABLE POSITION.

Area	How to Position the Straps		Straps to use
Torso (around waist)			
Torso (over shoulder)			
Hands & feet			
Arms & legs			
Face & exercise			

Hold, Hang or Lean
(no hanger provided)

(EN) Use & Care Guide

FOR YOUR SAFETY AND CONTINUED ENJOYMENT OF THIS PRODUCT, ALWAYS READ INSTRUCTION LEAFLET CAREFULLY BEFORE USING.

The product is intended for private use only and is not intended for commercial use or for use in other applications. This is not a medical product or a medical device.

SPECIFICATIONS:

Name: STYLPRO™ Red light Full Body panel

Model: RBP01A

Device Size: 400x 285 x 12.4 mm

Materials: Silicone

Input Voltage: 110V -240V

Output: 12v --- 2A

Weight: 565g

INCLUDES:

1 X STYLPRO™ Red Light Full Body Panel

2 X Adjustable body straps


1 X UK power supply

1 X EU plug adapter

1 X Protective eye glasses

1 X Glasses pouch

1 X Travel storage bag

12V  3A

BEFORE USE – IMPORTANT

LIGHT SENSITIVITY TEST (REQUIRED BEFORE FIRST USE)

BEFORE USING THIS DEVICE ON YOUR FACE OR BODY, PERFORM A LIGHT SENSITIVITY TEST:

1. PRESS THE ON/OFF BUTTON TO TURN THE DEVICE ON.
2. PLACE THE PANEL ON YOUR ARM AND RUN A 5-MINUTE SESSION.
3. AFTER THE SESSION, TURN THE DEVICE OFF AND CHECK THE AREA.

NOTE:

MILD REDNESS MAY APPEAR IMMEDIATELY AFTER USE DUE TO WARMTH FROM THE DEVICE AND SHOULD FADE QUICKLY.

IF REDNESS LASTS LONGER THAN 2 HOURS, DISCONTINUE USE. YOU MAY HAVE SENSITIVITY TO LIGHT - SEEK MEDICAL ADVICE BEFORE CONTINUING USE.

SAFETY INSTRUCTIONS:

1. THIS APPLIANCE IS DESIGNED FOR INDOOR HOME USE ONLY.

2. CONSULT A DOCTOR BEFORE USE, IF YOU HAVE ANY OF THE CONDITIONS BELOW, USE OF THE BODY & FACE SMOOTHER IS NOT RECOMMENDED WITHOUT YOUR DOCTOR'S APPROVAL: METAL IMPLANTS, RECENT SURGERY, OR HAVE A MEDICAL IMPLANT SUCH AS A PACEMAKER, CURRENTLY PREGNANT OR BREASTFEEDING, NAUSEA, HEADACHE, FEVER, SKIN DISEASES, BURNS OR REDDENING, BROKEN BONES, CANCER OR OTHER MALIGNANT CONDITIONS, OSTEOPOROSIS AND MENTAL ILLNESS.

3. FOR ADULT USE ONLY. DO NOT USE THE RED LIGHT BODY PANEL WITH AN INFANT, CHILD, PET, INCAPACITATED OR IMMOBILE PERSON, OR PERSONS WITH REDUCED PHYSICAL OR MENTAL CAPABILITIES WITHOUT SUPERVISION. THIS IS NOT A TOY.

4. IF THE RED LIGHT BODY PANEL MALFUNCTIONS OR IS DAMAGED IN ANY WAY, STOP USING IT IMMEDIATELY. CONTACT STYLPRO™ CUSTOMER SERVICE (SEE SECTION LIMITED TWO-YEAR WARRANTY FOR CONTACT DETAILS).

5. USE LED LIGHTS WITH CARE. EXTENDED USE AND PROLONGED EXPOSURE OF LED LIGHTS INTO THE EYES MAY CAUSE DAMAGE. LEDS ARE NOT REPLACEABLE.

6. DO NOT USE ON BROKEN, IRRITATED, SWOLLEN, SUNBURNT, OR NUMB/SENSITIVE SKIN, AREAS OF SWELLING, BRUISING, INFLAMMATION, OR SKIN CONDITIONS.

7. DO NOT LOOK DIRECTLY INTO THE LED LIGHTS.

8. AVOID USING WITH SERUMS CONTAINING LIGHT SENSITIVE INGREDIENTS OR POTENTIAL PHOTOSENSITIZERS SUCH AS VITAMIN C, AHAS, BHAS, ESSENTIAL OILS, BENZOYL PEROXIDE, AND RETINOLIDS.

TEST YOURSELF FOR LIGHT SENSITIVITY ON YOUR ARM PRIOR TO USING THE RED LIGHT BODY PANEL

HOW TO USE:

1. CONNECT THE DEVICE TO THE POWER CABLE, THEN PLUG THE POWER ADAPTER INTO A WALL SOCKET.
2. ENSURE THE AREA IS CLEAN, DRY AND FREE FROM PRODUCTS.

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
3. ALWAYS WEAR THE SUPPLIED EYE PROTECTION WHEN THE PANEL IS FACING OR POSITIONED NEAR THE FACE.
4. PREPARE THE PANEL FOR YOUR DESIRED SESSION POSITION:
 - THE PANEL CAN BE HELD, LEANED AGAINST A STABLE SURFACE, OR MOUNTED TO DIRECT LIGHT TOWARDS THE FACE OR BODY.
 - THE PANEL CAN ALSO BE WRAPPED AROUND THE BODY USING ONE OR TWO STRAPS, DEPENDING ON THE AREA BEING TREATED.
5. LONG PRESS THE ON/OFF BUTTON TO START THE SESSION.
 - FOUR INDICATOR LIGHTS WILL ILLUMINATE, EACH REPRESENTING 5 MINUTES.
 - ONE LIGHT WILL TURN OFF EVERY 5 MINUTES.
 - THE DEVICE WILL AUTOMATICALLY TURN OFF AFTER 20 MINUTES.
6. ENJOY YOUR SESSION.
WE RECOMMEND DAILY USE FOR 10-20 MINUTES PER BODY AREA AS PART OF YOUR REGULAR ROUTINE.

HOW TO WEAR:

STRAP SETUP (BEFORE POSITIONING)
ATTACH THE LONGER AND SHORTER STRAPS TO EACH SIDE OF THE PANEL, ENSURING THE SOFT (FLUFFY) VELCRO FACES OUTWARD.

TORSO (AROUND WAIST)

HOLD THE PANEL BY EACH STRAP AND POSITION IT AROUND YOUR WAIST.
FEED THE LONGER STRAP THROUGH THE PLASTIC LOOP OF THE SHORTER STRAP.
SECURE THE STRAP SO THE PANEL SITS COMFORTABLY AND LOOSELY AROUND THE BODY.

TORSO (OVER SHOULDER)

HOLD THE PANEL BY THE SHORTER STRAP AND LIFT IT OVER YOUR SHOULDER.
FEED THE LONGER STRAP THROUGH THE PLASTIC LOOP OF THE SHORTER STRAP.
SECURE SO THE PANEL RESTS COMFORTABLY ACROSS THE BODY.

HANDS & FEET

FEED THE LONGER STRAP THROUGH THE PLASTIC LOOP OF THE SHORTER STRAP AND TIGHTEN FULLY.
THIS FORMS A TUNNEL SHAPE, ALLOWING HANDS OR FEET TO SIT COMFORTABLY BENEATH THE PANEL.

ARMS & LEGS

REMOVE THE SHORTER STRAP.
LOOP THE LONGER STRAP THROUGH THE OPPOSITE SIDE OF THE PANEL AND TIGHTEN.
THE REMAINING STRAP LENGTH CAN BE WRAPPED AROUND THE PANEL TO SECURE IT IN PLACE.

EXERCISE OR FACE USE

REMOVE BOTH STRAPS.
HOLD THE PANEL IN YOUR HANDS, LEAN IT AGAINST A STABLE SURFACE, OR HANG IT USING A HANGER (NOT INCLUDED).
ALWAYS WEAR EYE PROTECTION IF THE PANEL IS FACING THE FACE.

HOW TO CLEAN:

1. USE A DISINFECTANT WIPE TO WIPE DOWN THE LED DEVICE. WE RECOMMEND TO CLEAN AFTER EVERY USE.

SAFETY CONSIDERATIONS:

Use of a Red Light skincare device in any way other than those described in this manual voids the warranty and could result in serious injury.

The STYLPRO™ THE RED LIGHT BODY PANEL is an electrical device. To avoid electric shock, adhere to the following safety considerations:

- Do not get your device wet or submerge your device in water.
- Do not wash electrical parts with water or liquids.
- Do not directly touch the device or switches with wet hands.
- Do not apply or touch the device to wet or damp skin.
- Do not use the device if it has been dropped, left outside, or has encountered water or liquids.
- Do not modify, or attempt to modify or repair, this equipment. There are no user serviceable parts.
- Do not use or store your device in hot temperatures that are or can exceed 54°C (130°F).

WARNING:

TO REDUCE THE RISK OF BURNS, ELECTROCUTION, FIRE, OR INJURY TO PERSONS:

1. READ ALL INSTRUCTIONS BEFORE USING THE APPLIANCE.
2. IF YOU FEEL ABNORMALITIES IN YOUR SKIN DURING USE OR AFTER USE, PLEASE CONSULT YOUR MEDICAL SPECIALIST.
3. TO PROTECT AGAINST ELECTRICAL HAZARDS, DO NOT USE WHILE BATHING OR IN THE SHOWER.
4. DO NOT OPERATE THIS APPLIANCE IF THE APPLIANCE MALFUNCTIONS OR HAS BEEN DROPPED OR DAMAGED IN ANY MANNER. RETURN APPLIANCE TO NEAREST SERVICE FACILITY FOR EXAMINATION, ELECTRICAL OR MECHANICAL ADJUSTMENT, OR REPAIR. NO REPAIRS SHOULD BE ATTEMPTED BY THE CONSUMER.
5. USE APPLIANCE ONLY FOR INTENDED USE AS DESCRIBED IN THIS INSTRUCTION BOOK. THE USE OF ATTACHMENTS NOT RECOMMENDED BY STYLPRO™ MAY CAUSE HAZARDS.
6. DO NOT PLACE ANY PART OF THE PRODUCT NEAR HEAT SOURCES SUCH AS RADIATORS, COOKERS, ETC DO NOT LAY ON A SOFT

- SURFACE SUCH AS A BED AS DEVICE MAY OVERHEAT.
- 7. KEEP THE DEVICE AWAY FROM CORROSIVE LIQUIDS.**
- 8. DO NOT SHARE THIS PRODUCT WITH ANYONE ELSE.**
- 9. DO NOT PLACE OR STORE WHERE IT CAN FALL OR BE PULLED INTO A TUB OR SINK. IF IT FALLS INTO WATER, UNPLUG IMMEDIATELY—DO NOT REACH INTO WATER.**
- 10. DO NOT TOUCH PLUG/CONTROLLER WITH WET HANDS OR USE ON WET/DAMP SKIN OR CLOTHING.**
- 11. DO NOT BUNCH, CREASE, OR TIGHTLY FOLD DURING OPERATION; DO NOT SIT/KNEEL ON IT OR PLACE HEAVY/SHARP OBJECTS ON IT.**
- 12. KEEP AWAY FROM AEROSOLS/SOLVENTS, OXYGEN USE, BEDDING, CURTAINS, OR OTHER FLAMMABLE MATERIALS.**
- 13. ALWAYS WEAR SUPPLIED EYE PROTECTION.**

ELECTRICAL SAFETY INFORMATION:

- 1.** Ensure the appliance is connected to a correctly rated power supply matching the voltage indicated on the product label.
- 2.** UK/EU-non UK customers: USE ONLY WITH A 230V(UK)/220-240-240V (EU-non UK) AC OUTLET. Use of the blanket in other countries with third party electrical transformers is at your own risk.
- 3.** Do not use with an extension cord. Do not use with generators, power converters or inverters. Use of the blanket in other countries with third party electrical transformers is at your own risk.
- 4.** Do not use if there is fraying of the wires, the AC cord or if the controller, cord or connectors appear mechanically or electrically compromised/damaged. If malfunction occurs, discontinue use immediately and contact STYLPRO™ Customer Service.
- 5.** The plug must be fitted with a fuse of the same rating (ASTA approved to BS 1362). If in doubt, consult a qualified electrician.
- 6.** Keep power cord away from heated surfaces; never pull or carry the device by its cord.
- 7.** Inspect plug, cord, and controller regularly. Do not use if damaged, overheated, dropped, or malfunctioning. No user-serviceable parts—do not open or modify.

STORAGE:

ENSURE DEVICE IS 'OFF' BEFORE STORAGE. ENSURE THE DEVICE IS UNPLUGGED FROM ANY POWER SUPPLY. PLACE IN COOL, DRY AREA OUT OF REACH OF CHILDREN, WHERE IT WILL NOT BE DAMAGED.

ALWAYS STORE THE DEVICE IN ITS NATURAL CURVED POSITION. ALLOW THE DEVICE TO COOL BEFORE STORING.

MAINTENANCE:

YOUR RED LIGHT BODY PANEL IS VIRTUALLY MAINTENANCE FREE. TO CLEAN TURN THE APPLIANCE OFF, LET THE DEVICE COOL TO ROOM TEMPERATURE, AND WIPE CLEAN WITH A DAMP CLOTH. IF ANY ABNORMAL CONDITION OCCURS, TURN OFF & CONTACT STYLPRO™ CUSTOMER SERVICES ON +44 (0) 20 3225 5559.

LIMITED TWO-YEAR WARRANTY:

OUR GOODS COME WITH GUARANTEES THAT CANNOT BE EXCLUDED UNDER THE EUROPEAN CONSUMER LAW. EU DIRECTIVE 1999/44/EC, THIS DIRECTIVE STATES 'A TWO-YEAR GUARANTEE APPLIES FOR THE SALE OF ALL CONSUMER GOODS EVERYWHERE IN THE EU. YOU ARE ENTITLED TO A REPLACEMENT OR REFUND FOR A MAJOR FAILURE YOU ARE ALSO ENTITLED TO HAVE THE GOODS REPAIRED OR REPLACED IF THE GOODS FAIL TO BE OF ACCEPTABLE QUALITY AND THE FAILURE DOES NOT AMOUNT TO A MAJOR FAILURE. STYLPRO™ WILL REPAIR OR REPLACE (AT OUR OPTION) YOUR APPLIANCE FREE OF CHARGE FOR TWENTY-FOUR MONTHS FROM THE DATE OF PURCHASE, IF THE APPLIANCE IS DEFECTIVE IN WORKMANSHIP OR MATERIALS. THIS WARRANTY IS AN ADDITIONAL BENEFIT AND DOES NOT AFFECT YOUR LEGAL RIGHTS. TO OBTAIN SERVICE UNDER THIS WARRANTY, PLEASE CONTACT STYLPRO™ CUSTOMER SERVICE AT TEAMTOM@STYLIDEAS.COM +44 (0) 20 3225 5559. TO MAKE A WARRANTY CLAIM YOU MUST RETAIN YOUR PROOF OF PURCHASE.

IMPORTANT: THIS WARRANTY DOES NOT COVER DAMAGE RESULTING FROM ACCIDENTS, MISUSE, ALTERATIONS, ABUSE, LACK OF REASONABLE CARE OR USE OTHER THAN DESCRIBED IN THESE INSTRUCTIONS, OR PERFORMANCE OF SERVICE OR REPAIRS BY AN UNAUTHORISED SERVICE CENTRE. THIS WARRANTY IS VALID ONLY ON MERCHANDISE SOLD BY AUTHORISED RETAILERS. NOTE, THE ITEM MUST BE THOROUGHLY CLEANED PRIOR TO RETURN.

CORRECT DISPOSAL OF THIS PRODUCT:

Disposal of Product: The WEEE symbol on the STYLIDEAS, the Original Accessories or packaging indicates that this device must not be treated as unsorted municipal waste but must be collected separately! Dispose of the STYLIDEAS Product and its Original Accessories via a collection point for the recycling of waste electrical and electronic equipment if you live within the EU and in other European countries that operate separate collection systems for waste electrical and electronic equipment. By disposing of the STYLIDEAS Product and its Original Accessories in the proper manner, you help to avoid possible hazards for the environment and public health that could otherwise be caused by improper treatment of waste equipment. The recycling of materials contributes to the conservation of natural resources. Therefore, do not dispose of your old electrical and electronic equipment with the unsorted municipal waste.

Disposal of Packaging Waste:

The packaging is made of environmentally friendly materials, which may be disposed through your local recycling facilities. By disposing of the packaging in the proper manner, you help to avoid possible hazards for the environment and public health. The symbol on the carton packaging indicates that the packaging is made of PAP.

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